


autoRM: an Automatic Retinotopic Mapping Tool for Mice

version 1.0

release date 2021/4/19

written by Ruix. Li

autoRM provides a fully automatic tool for mice retinotopic mapping. It can help you to locate the primary and several higher-order visual cortex of mice with calcium or optical imaging (figure from ref [1]).



What you need to use autoRM

1. **Psychopy3** (>=2020.1.3) to present visual stimulus
2. **NI-DAQ** USB-6501 digital I/O Device and NI-DAQmx driver for synchronization
3. **MATLAB** (>=2019a) with image processing toolbox


The content of autoRM

- **RetinotopicMappingv5.py** is a Python script to present visual stimulus for retinotopic mapping, a typical experiment last about 30 minutes.
- **RMDegMap.m** is a MATLAB function to calculate visual degree, naming **azimuth** and **elevation**.
- **RMSetParam.m** is a MATLAB app helps you to determine parameters used in **RMAreaMap.m**.
- **RMAreaMap.m** is a MATLAB function to identify visual areas, see ref [2] for details.
- **autoRM.m** is **the MATLAB function you use**, it calls **RMDegMap** **RMSetParam** and **RMAreaMap**, usually you don't need to use other functions.

How to use

1. MATLAB image processing toolbox is required.

2. Run `RetinotopicMappingv5.py` with psychopy3 before start recording. Recommend using a NI-DAQ digital I/O device to synchronize your camera with the visual stimulus.
3. The Python program will return the frame number required for recording. **Set the frame number in your camera control interface.**




4. After stimulus finish, `RetinotopicMappingv5.py` will save a txt log file and a json configuration file.

The json file is required in the following steps.


5. Convert your recording data into `.mat` data.
6. Add `RMDegMap.m`, `RMSetParam.m` and `RMAreaMap.m` to your MATLAB path, then call `autoRM` in MATLAB command window and follow the instructions to select recording data `.mat` file and `.json` configuration files.
7. A GUI will pop-up for you to adjust parameters.

If you use `autoRM` function, data will be automatically loaded and saved



Adjust parameters by looking at the figure. Close the window when finished.

8. `autoRM` will create a figure for you at the path of `.json` config file. This figure can be used as retinotopic mapping reference for your subsequential experiments.



Inputs and outputs

autoRM need the following input to work

- **XYT image data**: the image data for analysis
- **configuration json file**: the json file generated by `RetinotopicMappingv5.py`, which contains experiment configurations for data process.

autoRM save the following data

degMap mat file contains the following data

- **FOV**: a image of field of view
- **dataL2R, dataR2L**: trial average data of left to right and right to left
- **dataD2U, dataU2D**: trial average data of down to up and up to down
- **phaseMaps**: phase maps of trial average data obtained from FFT.
- **degMaps**: degree maps connect visual areas with visual field.
- **degMapAzi**: azimuth map. Here assume the azimuth at front of mouse is 0. Azimuth usually ranges from 0 to 120 degree.
- **degMapElv**: elevation map. Here assume the elevation at front of mouse is 0. Elevation ususally ranges from -40 to 40 degree.

areaMap mat file contains the following data

- **rSignMap**: raw sign map calculated from **degMapAzi** and **degMapElv**
- **signMap**: sign map after image process to reduce noise and artifacts.
- **areaMap**: a bit map indicate visal areas.
- **signFOV**: align sign map to FOV.
- **hltFOV**: high light visual areas in FOV.

autoRM save 2 figures

Endnote

References:

- [1] Marshel, James H., et al. "Functional specialization of seven mouse visual cortical areas." Neuron 72.6 (2011): 1040-1054.
- [2] Juavinett, Ashley L., et al. "Automated identification of mouse visual areas with intrinsic signal imaging." Nature protocols 12.1 (2017): 32.
- [3] Zhuang, Jun, et al. "An extended retinotopic map of mouse cortex." Elife 6 (2017): e18372.